VERMEIREN

Eclips X4 KIDS

INSTRUCTION MANUAL
MODE D'EMPLOI
GEBRUIKSAANWIJZING
GEBRAUCHSANWEISUN
GISTRUZIONI PER L'USO
MANUAL DE INSTRUCCIONES
INSTRUKCJA OBSŁUGI
NÁVOD K OBSLUZE

EN

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Contents

Preface	2
1 Product description	3
1.1 Intended Use.....	3
1.2 Technical specifications	4
1.3 Drawing	5
1.4 Accessories.....	5
1.5 Location identification plate.....	5
1.7 Explanation of symbols	6
1.8 Safety rules	6
2 Use	7
2.1 Carrying the wheelchair	7
2.2 Mounting the rear wheels (if quick release axle is available).....	7
2.3 Unfolding the wheelchair	7
2.4 Mounting or removing of the foot rests	7
2.5 Operating the brakes.....	8
2.6 Transfer in and out the wheelchair	8
2.7 Correct position in the wheelchair	8
2.8 Riding the wheelchair.....	9
2.9 Moving on slopes	9
2.10 Negotiating steps or curbs	9
2.11 Fold up the wheelchair	11
2.12 Taking off the wheels	11
2.13 Push bar	11
2.14 Transport in the car.....	12
2.15 Use of the wheelchair as seat in a motor vehicle (Eclips X4 kids)	12
3 Installation and adjustment	14
3.1 Tools	14
3.2 Manner of delivery	14
3.3 Adjusting the seat height and seat angle	15
3.4 Flexible seat.....	17
3.5 Anatomical seat.....	17
3.6 Adjusting seat depth	18
3.7 Flexible backrest	18
3.8 Anatomical backrest.....	19
3.9 Adjust the brakes	19
3.10 Adjusting of the foot rests	19
3.11 Adjusting of the armpads	20
3.12 Adjusting of the tip cap	21
3.13 Adjusting the stability and the maneuverability.....	21
4 Maintenance	22

Preface

First of all we want to thank you for putting your trust in us by selecting one of our wheelchairs.

The Vermeiren wheelchairs are the result of many years of research and experience. During the development, special attention was given to the ease of use and the serviceability of the wheelchair.

The expected lifetime of your wheelchair is strongly influenced by the care and maintenance of the wheelchair.

This manual will help you get acquainted with the operation of your wheelchair.

Following the user instructions and the maintenance instructions are an essential part of the warranty.

This manual reflects the latest product developments. Vermeiren has the right to introduce changes without the obligation to adapt or replace previously delivered models.

For any further questions, please consult your specialist dealer.

1 Product description

1.1 *Intended Use*

The wheelchair is intended for people with walking difficulties or no walking abilities.

The wheelchair is designed to transport 1 person.

The wheelchair is suited for indoor and outdoor use.

The user can propel the wheelchair by himself or have the wheelchair pushed by an attendant.

The different types of fittings and accessories, and the modular construction allow full use by persons disabled by:

- paralysis
- loss of limbs (leg amputation)
- limb defects or deformations
- stiff or damaged joints
- heart insufficiencies and poor blood circulation
- balance disturbances
- cachexia (decrease in muscle)
- and also for aged persons.

When providing for individual requirements:

- body size and weight (max. 80 kg)
- physical and psychological condition
- residential circumstances
- environment

should be taken into consideration.

Your wheelchair should only be used on surfaces where all four wheels are touching the ground and where there is sufficient contact to propel the wheels equally.

You should practice for use on uneven surfaces (cobblestones, etc.), slopes, curves and to get past obstacles (curbs, etc.).

The wheelchair should not be used as a ladder, nor is it a transport for heavy or hot objects.

When used on mats, carpeted floors or loose floor coverings, the floor covering can get damaged.

Use only Vermeiren approved accessories.

The manufacturer is not liable for damage caused by the lack of or improper service or as a result of not following instructions from this manual.

Compliance with the user and maintenance instructions are an essential part of the guarantee conditions.

1.2 Technical specifications

Technical terms below are valid for the wheelchair in standard settings.

If other foot rests / arm rests or other accessories are used, the tabulated values will change.

Make	Vermeiren	
Address	Vermeirenplein 1/15, B-2920 Kalmthout	
Type	Manual wheelchair	
Model	Eclips X4 kids	
Maximum occupant mass	80 kg	
Description	Dimensions	
Effective seat width	320 mm	370 mm
Overall width (depends on the seat width)	500 mm	550 mm
Effective seat depth	310 mm – 370 mm	330 mm - 390 mm
Description	Dimensions Min.	Dimensions Max.
Overall length with foot rest	830 mm	1000 mm
Folded length without foot rest	630 mm	800 mm
Folded width	310 mm	
Folded height	820 mm	980 mm
Total mass	± 15,00 kg	
Mass of heaviest part	9,40 kg	
Masses of parts that can be dismantled or removed	Foot rests: 1,80 kg; Rear wheels: 3,80 kg	
Static stability downhill	10° (in standard configuration)	
Static stability uphill	10° (in standard configuration)	
Static stability sideways	>16° (in standard configuration)	
Obstacle climbing	60 mm	
Seat plane angle	5°	
Seat surface height at front edge	420 mm	470 mm
Backrest angle	5°	
Backrest height	350 mm	400 mm
Distance between foot rest and seat	310 mm	
Angle between seat and foot rest	11°	
Distance between armpad and seat	150 mm	180 mm
Front location of arm rest structure	150 mm	
Handrim diameter	480 mm	
Horizontal location of axle (deflection)	70 mm	120 mm
Minimum turning radius	1530 mm	
Diameter PU Rear wheels	22"	
Tyre pressure, rear (driving) wheels (only for air tyres)	Max. 3.8 bar	
Diameter PU steering wheels	150 mm	
Tyre pressure, steering wheels (only for air tyres)	Max. 2.5 bar	
Storage and use temperature	+ 5 °C	+ 41 °C
Storage and use humidity	30%	70%

We reserve the right to introduce technical changes. Measurement tolerance ± 15 mm / 1,5 kg / °

Table 1: Technical specifications Eclips X4 kids


The wheelchair complies to the requirements set up in:

ISO 7176-8: Requirements and test methods for static, impact and fatigue strengths.

ISO 7176-16: Resistance to ignition of upholstered parts

ISO 7176-19: Wheeled mobility devices for use as seats in motor vehicles.

1.3 Drawing


1.4 Accessories

The following accessories are available for the Eclips X4 kids:

- Anterior pelvic belt (B58) for mounting on the tubes of the back
- Anti-tipping device (B78) for mounting on the bottom frame
- Arm rests (B05, B82)
- Foot rests (B06, BA6, BZ8, BZ7, BK6)
- Table (B12)
- Transit wheelchair (T30)

Contact your specialist dealer about other accessories. He will gladly advise you.

1.5 Location identification plate

1.7 Explanation of symbols

Maximum mass

Indoor and outdoor use

Maximum Safe Slope

CE conformity

Type designation

Not intended to be used as a seat in a Motor vehicle

1.8 Safety rules

- ⚠ To prevent injury and/or damage to your wheelchair, make sure that no objects and/or body parts are caught in the spokes of the driving wheels.
- ⚠ The parking brakes should be applied before getting into and out of the wheelchair.
- ⚠ When getting into and out of the wheelchair, do not stand on the footplates. These should be folded up beforehand, or the foot rests should be swung completely outwards out of the way.
- ⚠ Investigate the effects of shifting the centre of gravity on the behavior of the wheelchair, for example on up or down gradients, on laterally sloping ground, or when overcoming obstacles. Obtain support from an attendant.
- ⚠ If you want to pick up something (lying in front of, on the side, or to the rear of the wheelchair), you should not lean too far out to avoid tipping over.
- ⚠ When moving through doors, arches, etc. ensure that there is enough room at the sides so that you do not get your hands or arms caught or crushed and that there is no damage to the wheelchair.
- ⚠ Only use your wheelchair according to regulations. For example, avoid uncontrolled rolling against obstacles (steps, curbs, doorframes, etc) or dropping down from ledges. The manufacturer cannot assume liability for damage caused by overloading, collision or other improper use.
- ⚠ Stairs may only be negotiated when aided by another person. If furnishings like drive-on, ramps or lifts are available, use them.
- ⚠ When moving on the public highway, you are subject to the rules of the road.
- ⚠ When driving your wheelchair, you should not be under the influence of alcohol or medicine as in the case of driving other vehicles. This also applies to indoor driving.
- ⚠ When travelling outdoors, adapt your driving to weather and traffic conditions.
- ⚠ When transporting the wheelchair, never pick it up by grasping movable parts (arm rests, foot rests, etc.).
- ⚠ To be better visible when driving in the dark, wear the brightest possible clothing or clothes with reflectors, and check that the reflectors mounted on the sides and rear of the wheelchair are clearly visible.
- ⚠ Never exceed the maximum load of 80 kg.

2 Use

This chapter describes the everyday use. **These instructions are for the user and the specialist dealer.**


The wheelchair is delivered fully assembled by your specialist dealer. The instructions intended for the specialist dealer on how to set up the wheelchair are given in § 3.

2.1 Carrying the wheelchair

The best way to carry the wheelchair is to make use of the wheels and roll the wheelchair.

If this is not possible (e.g. when the rear wheels are taken off for transportation in a car), firmly grasp the frame on the front and the handgrips. Do not use the foot rests, arm rests or the wheels to grasp the wheelchair.

2.2 Mounting the rear wheels (if quick release axle is available)


1. Take the rear wheel and push on button ①.
2. Keep the button pushed in and mount the rear wheels axle till it stops.
3. Release the button ①.
4. Check that the wheel is secured.

2.3 Unfolding the wheelchair

⚠ CAUTION: Risk of clamping – Keep fingers away from moving parts of the wheelchair.

1. Position yourself behind the wheelchair.
2. Use the handgrips to open the wheelchair as much as possible.
3. Position yourself at the front of the wheelchair.
4. Push both seating tubes down till they are fixed in their position.
5. When the push bar is available mount the push bar to the right handgrip.

2.4 Mounting or removing of the foot rests

The mounting of the foot rests is done as follows:


1. Hold the foot rest sideways at the outside of the wheelchair frame and mount the tube hood ① into the frame.
2. Swing the foot rest inwards till it clicks in position.
3. Swing the footplate downwards.

To take off the foot rests:

1. Pull lever ②.
2. Swing the foot rest to the outside of the wheelchair till it comes loose from the guidance.
3. Pull the foot rest tube hood ① out the wheelchair frame.

2.5 Operating the brakes

- ⚠ **WARNING:** The brakes are not used to slow down the wheelchair during movements – Use the brakes only to prevent the wheelchair from unintended movements.
- ⚠ **WARNING:** Good operation of the brakes is influenced by wear and contamination of the tires (water, oil, mud, ...) – Check the condition of the tires before each use.
- ⚠ **WARNING:** The brakes are adjustable and can wear – Check the operation of the brakes before each use.

To apply the brakes:

1. Push the brake levers ① forward till you feel a distinctive click.


- ⚠ **CAUTION:** Risk of unintended movement – Make sure the wheelchair is on a flat horizontal surface before releasing the brakes. Never release both brakes simultaneously.

To release the brakes:

1. Release one brake by pulling the lever ① backwards.
2. Hold the hand-rim of the released wheel with your hand.
3. Release the second brake by pulling the lever backwards.

2.6 Transfer in and out the wheelchair

- ⚠ **CAUTION:** In case you cannot perform the transfer in a safe manner, ask someone to assist you.
- ⚠ **CAUTION:** Risk of tipping over of the wheelchair – Do not stand on the footplates.

1. Position the wheelchair as close as possible to the chair, couch or bed to/from you wish to transfer.
2. Check both brakes from the wheelchair are in the on position.
3. Fold the footplates upwards to prevent standing on them.
4. Transfer to/from the wheelchair or let your parent or attendant lift you out the wheelchair.

2.7 Correct position in the wheelchair

Some recommendations for a comfortable use of the wheelchair:

- Position your backside as close as possible to the backrest.
- Make sure your upper legs are horizontal – If needed adjust the length of the foot rests. (see § 3.10.1)

2.8 Riding the wheelchair

- ⚠ **WARNING:** Risk of clamping – Prevent your fingers from being caught by the wheels spokes.
 - ⚠ **WARNING:** Risk of clamping – Be careful passing through restricted passages (e.g. doors).
 - ⚠ **WARNING:** Risk of burns – Be careful when driving in hot or cold environments (sunshine, extreme cold, saunas, etc.) for a long amount of time. Touching surfaces which can assume the environment temperatures, can cause burns.
1. Release the brakes.
 2. Take both hand rims at their highest position.
 3. Lean forward and push the hand rims forward until straight arms.
 4. Swing your arms loosely back to the hand rims topside and repeat the movement.

2.9 Moving on slopes

- ⚠ **WARNING:** Control your speed – Moving on slopes as slow as possible.
- ⚠ **WARNING:** Consider the capacities of your attendant – If your attendant does not have enough force to control the wheelchair, put on the brakes.
- ⚠ **WARNING:** Risk of tipping over – Lean forward to move your centre of gravity forward. To improve a better stability.

1. If available on the wheelchair, wear the safety belt.
2. Do not attempt moving on too high slopes. The maximum slope angles (upwards and downwards) are mentioned in table 1 for Eclips X4 kids.
3. Ask an attendant to help you moving on the slope.
4. Lean forward to move your centre of gravity forward.

2.10 Negotiating steps or curbs

2.10.1 Getting down steps or curbs

Riding down low curbs can be done moving forwards. Make sure that the footplates do not touch the ground.

A practiced user can negotiate small steps or curbs by himself:

- ⚠ **WARNING:** Risk of tipping over – If you do not have enough experience with your wheelchair, ask assistance of an attendant.
1. Bring balance on the rear wheels to reduce the pressure on the front wheels.
 2. Negotiate the curbs.

Higher curbs can be taken forward with an attendant:

1. Ask the attendant to tip the wheelchair slightly backwards.
2. Get past the curbs while moving on the rear wheels.
3. Put the wheelchair back on the four wheels.

An experienced user can negotiate higher curbs by himself. This is best done backwards.

1. Turn the wheelchair to have the rear wheels facing the curb.
2. Lean forward to move your centre of gravity forward.
3. Move the wheelchair close to the curbs.
4. Use the hand rim to roll-off wheelchair from the curb in a controlled manner.

2.10.2 Moving up steps or curbs

Moving up steps or curbs with attendant as follows:

1. Prevent the footplates from touching the curb.
2. Ask the attendant to tip the wheelchair backwards, just enough to move the front wheels over the curb.
3. Lean backwards to move your centre of gravity above the rear wheel.
4. Place the front wheels on the curb.
5. Roll rear wheels of the wheelchair over the curb.

Higher curbs are negotiated backwards:

1. Turn the wheelchair to have the rear wheels facing the curb.
2. Lean backwards and move your centre of gravity above the rear wheels.
3. Ask the attendant to pull the wheelchair on the curb.
4. Take back the normal position in the wheelchair.

An experienced user can negotiate curbs by himself:

⚠ WARNING: Risk of tipping over – If you have not enough experience to control the wheelchair, get help from an attendant.

1. Drive until the curbs.
2. Ensure that the foot rest do not touch the curbs.
3. Lean backwards so you are balancing on the rear wheels.

4. Role the front wheels balancing over the curbs.
5. Bend forward for more stability.
6. Role the rear wheels over the curbs.

2.10.3 Taking of stairs

Taking of stairs while you staying in the wheelchair shall be according following rules:

⚠ WARNING: Risk of tipping over – Taking of stairs shall always with 2 attendants.

1. Remove the foot rests.
2. One attendant tip the wheelchair slightly backwards.
3. The second attendant take the front of the frame.
4. Stay calm, avoid sudden movements and keep your arms inside the wheelchair.
5. Take the steps on the rear wheels of the wheelchair.
6. Mounting the foot rests back after taking the stair.


2.11 Fold up the wheelchair

⚠ CAUTION: Chance of pinching – Do not place fingers between the components of the wheelchair.

1. When the push bar is available remove the push bar.
2. Fold or remove the footplates (see § 2.4).
3. Take the seat on the front side and backside and pull it up.
4. Push on the handgrips to fold up the wheelchair further.

2.12 Taking off the wheels

To facilitate the transport off the wheelchair the rear wheels can be taken off:

1. Make sure the brakes standing in the off position.
2. Take the wheelchair to the side frame where you want to remove the wheel.
3. Press the button ① in the center of the wheel hub.
4. Pull the wheel away from the frame.

2.13 Push bar


⚠ WARNING: Risk of injury - Make sure that all screws are properly tightened with the appropriate tool.

⚠ WARNING: Risk of injury - When equipped with the push bar, the push bar shall always be mounted (→ to increase the stability of the backrest).

Your wheelchair can have a push bar. The push bar increase the stability of the backrest and offer you the possibility of pushing the wheelchair by an attendant. To do this, use the pre-mounted handgrips.

As there are different push bars for different seat widths, take care that no tension between the handgrips of the back tubes is created.

The mounting of the push bar is done as follows:

1. Secure the push bar ① by fastening the left side to the end of back tube handgrip ② with the Philips screw ③ in the applicable sunk threads.
2. Screw gently the Philips screw ③ on the right hand side into the thread of the back tube.
3. Swing the push bar to the right.
4. Hook the suspension of the push bar ① over the screw thread.

Use the two holes on the left handgrip to adjust the tension of the backrest.

To take off the push bar:

1. Loosen the Philips screw ③ on the right handgrip.
2. Remove the hook from the right handgrip.

2.14 Transport in the car

- ⚠ **WARNING: Risk of injury – See that the wheelchair is attached properly. So you can avoid injury from the passengers during collision or sudden braking.**
- ⚠ **WARNING: Risk of injury – Use for attaching the wheelchair and passenger NEVER the same seatbelt.**

To transport the wheelchair in the car use following steps:

1. Remove foot rests and accessories.
2. Store foot rests and accessories safely.
3. If possible, fold the wheelchair and remove the wheels.
4. Place the wheelchair in the luggage place.
5. If the wheelchair and the passenger compartment is NOT separated, attach the frame of the wheelchair securely to the vehicle. You can use the available safety belts in the vehicle.

2.15 Use of the wheelchair as seat in a motor vehicle (Eclips X4 kids)

- ⚠ **WARNING: The wheelchair has passed the crash test of ISO 7176-19: 2008 and, as such, has been designed and tested for use only as forward-facing seat in a motor vehicle.**
- ⚠ **WARNING: The wheelchair's pelvic belt alone is not suited as an occupant restraint belt.**

The wheelchair is tested using the four-point strap-tie system and a 3-point occupant-restraint system.

Whenever feasible, use the seat of the vehicle and store the wheelchair in the cargo area.

Steps to secure the wheelchair in a motor vehicle:

1. Check that the vehicle is equipped with a suitable wheelchair tie down and occupant-restraint system, conform ISO 10542.
2. Check that the components of the wheelchair tie down and occupant restraint system are not frayed, contaminated, damaged or broken.
3. If equipped with an adjustable seat and/or back tilt, make sure that the wheelchair user is sitting as upright as possible. If the user's condition prevents this, a risk assessment should be done to evaluate the user's safety during transit.
4. Remove all mounted accessories such as trays and respiratory equipment, and secure them in a safe place.
5. Position the wheelchair facing forward in the travelling direction, centrally between the tie-down rails mounted in the floor of the vehicle.
6. Make sure that the indicated zones around the wheelchair user are clear from rigid vehicle parts.


figure 1


figure 2

7. Mount the front securement straps according to the instructions of the strap-system manufacturer at the indicated place. (figure 3)
This place is marked on the wheelchair with a symbol. (figure 4)
8. Roll back the wheelchair until the front straps are tight.
9. Apply the wheelchair brake.
10. Mount the back securement straps according to the instructions of the strap-system manufacturer at the indicated place. (figure 3)
11. This place is marked on the wheelchair with a symbol. (figure 4)



figure 3


figure 4

Steps to secure the wheelchair user:

1. Remove both arm rests.
2. If present, attach the wheelchair's pelvic belt.
3. Attach the occupant restraint belts according to the instructions of the strap-system manufacturer.
4. Wear the pelvic belt low across the front of the pelvis, so that the angle of the pelvic belt is within the preferred zone of 30° to 75° to the horizontal, similar to that shown below.

5. A steeper (greater) angle within the preferred zone is desirable.
6. Adjust the belt tightly according to the instructions of the strap-system manufacturer, consistent with the user's comfort.
7. Ensure that the restraint belt connects in a straight line to the anchor point in the vehicle and that no bends in the belt are visible, for instance at the axle of the rear wheel.
8. Install the arm rests, if desired. make sure that belts are not twisted or held away from the body by wheelchair components such as arm rests or wheels.

3 Installation and adjustment

The instructions in this chapter are for the specialist dealer.

The Vermeiren Lightweight wheelchairs Eclips X4 kids has been designed entirely for your comfort. This wheelchair can be adjusted to individual requirements as explained below.

To find a service facility or specialist dealer near you, contact the nearest Vermeiren facility. A list of Vermeiren facilities can be found on the last page.

- ⚠ **WARNING:** Risk of unsafe settings - Use only the settings described in this manual.
- ⚠ **WARNING:** Variation of allowed adjustments can still change the stability of your wheelchair (tilt back or sideways).

3.1 Tools

To set up the wheelchair the following tools are needed.

- Wrench set n° 10 to n° 19
- Allen key set n° 4 to n° 5
- Screwdriver Phillips head

3.2 Manner of delivery

The Vermeiren Eclips X4 kids shall be delivered with:

- 1 frame with arm rests (cloth protectors), rear and front wheels
- 1 pair foot rests
- Tools
- Manual
- Accessories


3.3 Adjusting the seat height and seat angle

- ⚠ **WARNING:** Risk of injury - Make sure that all screws are properly tightened with the appropriate tool.
- ⚠ **WARNING:** Risk of tipping over - When the driving wheels are adjusted, the position of the wheels are altered.

The Vermeiren Eclips X4 kids is adjustable in different seat heights and different seat angles by changing the wheels position. At each height there is a different setting of the front and rear wheels.

Change the seat height and seat angle according following steps:

* *Front wheel:*


1. Loosen the screws (A) from the steering wheel adapter (B).
 2. Choose correct size for the adapter (B) and fork.
 3. Place the steering wheel adapter in the desired position (2 different positions: steps from 34 mm).
 4. Tighten the screws (A) properly.
- or
5. Loosen the screw (D) from the front fork (E).
 6. Choose the correct position. **Do not use the holes (4) and (5).** (3 different positions: steps from 12 mm).
 7. Tighten the screws (D) properly.


Adjustment of the slope nuts (C) from the steering wheels.


- ⚠ **WARNING:** Risk of injury - Make sure that the eccentric screw settings are identical on both front-wheel adapters and, before using the wheelchair.

This figure illustrates the principle. The position of the screw thread is indicated by a number.


- A** Shows the frame when the seat is horizontal and the adapter is perpendicular to the ground.
- B** Shows the frame when the seat is tilted slightly and the adapter is not perpendicular.


** Rear wheel:*

1. Remove the rear wheels (see § 2.12).
2. Loosen the fixing screws (A) from the rear axle block (B).
3. Push the axle block (B) upwards or downwards to the desired height (steps from 15 mm).
4. Tighten the screws (A) properly.
5. Mount the rear wheels.
6. Adjust the brakes according § 3.9.

It is also possible to rotate the rear axle block (B) (2 positions: steps from 15 mm).

Some of the wheelchairs Eclips X4 kids Aktief are made with an adjustable axle adjuster. With these adjuster you can also change the seat height.


1. Remove the screws ① and ②.
2. Push/Pull the adjustable axle adjuster ③ upwards or downwards to the required height.
3. Retighten the screws ① and ②.

3.4 Flexible seat

⚠ WARNING: Risk of injury - Make sure that all screws are properly tightened with the appropriate tool.

The seat frame is available in different seat depths.

To close the gap between the back frame and the seat, fold the flap of the backrest under the seat (Velcro fastening).


Changing the seat cover as follows:

1. Remove the screw at the front end of the seat frame.
2. Pull the locking pin © out.
3. Now the seat © can be removed by pulling it gently forwards out the sleeve of tube ©.

Do the reverse when you want to replace the seat.

3.5 Anatomical seat


If you have opted for an anatomical padded seat, place it on the tubular cross-hinge so that the seat's supports engage equally with the cross-hinge. In the wooden base you will find slots that are placed in the hinge's cross-pieces. To fold up the wheelchair, remove the seat by pulling it upwards.

3.6 Adjusting seat depth

⚠ WARNING: Risk of injury - Make sure that all screws are properly tightened before using the wheelchair.

The seat depth can be adjusted from 310 mm to 370 mm in steps from 20 mm for seat width 320 mm and from 330 mm to 390 mm in steps from 20 mm for seat width 370 mm.


Adjust the seat depth according following procedure:

1. Remove the arm rests (cloth protectors), foot rests and rear wheels.
2. Loosen the screws ④.
3. Loosen the screws ⑤.
4. Remove the crosses from the seat located on the main frame ③ by folding the wheelchair.
5. Remove the screws ⑥ located in the blocks for the cross ⑦ on the main frame ③.

6. Remove the adjustable axle adapter, if available on your wheelchair.
7. Adjust the seat depth to the desired position by moving the back frame ① and/or seat frame ② steps from 20 mm.
8. Tighten all the screws properly by hand.
9. Fold the flap of the backrest under the seat (Velcro fastening), to close the gap between the back frame and the seat.
10. Adjusting the seat is not necessary.
11. Replace the arm rests (cloth protectors), foot rests and rear wheels. Mount the arm rests in the correct hole of attachment plate ⑧.

3.7 Flexible backrest

⚠ WARNING: Risk of injury - Make sure that all screws are properly tightened before using the wheelchair.


The height of the backrest can be adjusted by setting the entire back tube.

For adjusting the height of the back cover, loosen the fixing screws ① behind the back tube. Then you can move the cover upwards or downwards. There are fixing holes every 20 mm in the back tube where the screws can be retightened.

3.8 Anatomical backrest

- ⚠ **WARNING:** Risk of injury – Check that the backrest is secured before use.
- ⚠ **WARNING:** Risk of injury – Avoid getting your fingers caught during the assembly process.


To use the anatomical backrest, proceed as follows:

1. Unscrew the flexible backrest.
2. Loosely screw on the brackets for the backrest.
3. Insert the backrest into the brackets from above.
4. Tighten the screws holding the brackets for the backrest.
5. Insert and tighten the safety screws for the backrest.


3.9 Adjust the brakes

- ⚠ **WARNING:** Risk of injury – Brakes may only be adjusted by your specialist dealer.

Adjust the brakes according following rules:

1. Install the wheels according § 2.2.
2. Disconnect the brakes by pulling lever ① backwards.

3. Loosen the bolts ② so the brake mechanism can slide over the guide ③.
4. Pull the brake mechanism over the guide ③ to the desired position.
5. Retighten the bolts ②.
6. Check working of the brakes.
7. If necessary repeat the above steps until the brakes are adjusted well.

3.10 Adjusting of the foot rests

3.10.1 Length of the foot rests

- ⚠ **CAUTION:** Risk of damage – Avoid that the foot rests make any contact with the ground. Keep a minimum distance from 60 mm above the ground.


Adjust the length of the foot rests as follow (2 positions possible: steps 35 mm → holes outside tube, 4 positions possible: steps 25 mm → inner tube) (Fig A):

1. Remove the screw ①.
2. Adjust the length of the foot rest to a comfortable length.
3. Tighten the screw ① properly.

3.10.2 Depth of the footplate


Adjust the depth of the footplate as follow (in three different positions: steps 20 mm):

1. Remove the screws ① (Fig B).
2. Adjust the depth of the foot rest to a comfortable position (Fig A).
Position 1: hole 2 and 4
Position 2: hole 3 and 5
Position 3: hole 4 and 6
3. Tighten the screws ① properly (Fig B).

3.11 Adjusting of the armpads

⚠ CAUTION: Risk of clamping – Do not place fingers, buckles and clothes between the splash guards and rear wheels.

The armpads of the wheelchair can be adjusted in 2 different positions.

Height of armpad	Holes
150 mm	A
180 mm	B

Table 2: Height of armpad


1. Remove the screws ①.
2. Grasp the cloth protector ② on the top.
3. Remove the cloth protector ②.
4. Adjust the cloth protector ② in the desired position.
5. Mount the cloth protector ②.
6. Tighten the screws ① properly.
7. Check that the cloth protector ② is firmly fixated.

3.12 Adjusting of the tip cap

⚠ WARNING: Risk of injury - Make sure that all screws are properly tightened with the appropriate tool.

A tip cap is mounted on the lower back frame that can be adjusted in its height depending on the seat height of the wheelchair.

The tip cap can be placed in 2 different position with a difference of 15 mm.

1. Remove the screws ①.
2. Adjust the height of the tip cap to a comfortable position.

Tip cap height	Tube tip cap ①	Tube frame ②
115 mm (Standard)	Hole 1 and 3	Hole 1 and 3
130 mm	Hole 1 and 3	Hole 2 and 4

Table 3: Tip cap heights for wheelchair with seat height of 420 mm

3. Tighten the screws properly ①. (Fig C)


If an anti-tipping is installed, it replaces the addition of the tip cap and can be used accordingly.

3.13 Adjusting the stability and the maneuverability

⚠ WARNING: Risk of injury - Make sure that all screws are properly tightened with the appropriate tool.

⚠ WARNING: Risk of injury - Make sure that the axle adapters ① of the rear wheels on both sides are mounted in the same position.


⚠ WARNING: Risk of tipping over - Be aware that you change the stability of the wheelchair.

The driving wheels are mounted in an axle adapter ① at the bottom of the back frame ②.

1. Loosen the fixing screws ③.
2. Turn the axle adapter ① to the back until the axle mounting is behind the back frame.
3. Tighten the screws properly ③.
4. Adjust the brakes.

Some of the wheelchairs are made with an adjustable axle adjuster. With these adjuster you can also change the stability and the manoeuvrability of your wheelchair.

1. Remove the screws ① and ②.
2. Push/Pull the adjustable axle adjuster ③ upwards or downwards to the required height or mount the rear wheel in another axle mounting ③.
3. Retighten the screws ① and ②.
4. Adjust the brakes.

Please note that your wheelchair might become tipping over. To prevent this, you can also mount an anti-tipping device.

4 Maintenance

For the maintenance manual of the manual wheelchairs refer to the Vermeiren website: www.vermeiren.com.

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CONTRACTUAL WARRANTY

We offer 5 years of warranty on standard wheelchairs, lightweight wheelchairs 4 years. Electronic wheelchairs, tricycles, beds and other products: 2 years (batteries 6 months) and multipoosition wheelchairs 3 years. This warranty is limited to the replacement of defective or spare parts.

APPLICATION CONDITIONS

In order to claim its warranty, part "B" of this card has to be given to your official Vermeiren dealer. The warranty is only valid when parts are replaced by Vermeiren in Belgium.

EXCEPTIONS

This warranty is not valid in case of:

- damage due to incorrect usage of the wheelchair,
- damage during transport,
- involvement in an accident,
- a dismount, modification or repair carried outside of our company and/or official Vermeiren dealership,
- normal wear of the wheelchair,
- non-return of the warranty card

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N.V. VERMEIREN N.V.
Vermeirenplein 1/15
B-2920 Kalmarhout
BELGIUM

Notes

A series of horizontal dotted lines for taking notes.

Notes

A series of horizontal dotted lines for writing notes.

Vermeiren GROUP
Vermeirenplein 1 / 15
B-2920 Kalmthout
BE

website: www.vermeiren.com

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All rights reserved, including translation.