

708 Kids

Manual kids wheelchair

Compact sturdy wheelchair in steel. Offset of the rear wheels for improved stability. Equipped with the Quickrelease system for fast and easy removal of the rear wheels.



VERMEIREN
we care for you

Mobility

Our mobility can challenge us, even from a young age. At **VERMEIREN** we try to make sure that you can overcome these challenges together, no matter short- or long-term.

VERMEIREN offers a wide range of products for children and young adults to ensure the best mobility solution. From care wheelchairs to flashy active wheelchairs, our products offer something for everyone.

Colours



Robust push brakes



Tipping piece

Technical data

Total width	500 / 550 mm
Total height	905 mm
Total length	985 mm
Folded width	310 mm
Seat width	320 / 370 mm
Height armrests	210 mm

Products designed to suit you



Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **708 Kids** complies with the highest norms, laid down by certified test institutes in most European countries.



Seat height	445 mm
Height backrest	420 mm
Seat depth	350 mm
Weight	18 kg
Maximum user weight	60 kg
Additional info	 

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you