

Goliat Artro

Heavy-duty rollator

Heavy-duty upright rollator with forearm supports.

Goliat Artro has a modern and strong design with height-adjustable forearm supports to relieve stress on hands and wrists whilst standing upright. With a maximum user weight of 200 Kg and a wider frame, this rollator is an excellent choice for users who need extra stability and sturdiness when navigating. The forearm supports provide an upright position to relieve joints and improve posture.





Extra strong bolts for more safety and stability



A hidden shopping compartment (max. 5 Kg)



A strong and well-made foam seat



Detail crutch holder

Colour(s)



red



silver

Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, **Goliat Artro** complies with the highest norms, laid down by certified test institutes in most European countries.

Technical data

■ Total width	650 mm
■ Total height	1140 - 1310 mm
■ Total length	750 mm
■ Seat width	470 mm
■ Seat height	490 - 560 mm
■ Seat depth	410 mm

■ Multifunctional brakes with lock	
■ Wheel diameter	8"
■ Weight	12,35 kg
■ Maximum user weight	200 kg
■ Maximum load basket	5 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.



VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you