

28

Double cross

This chair with double cross and fixed back can support a weight up to 110 kg (150 kg for wheelchairs with anatomic cushions).



VERMEIREN
we care for you

Solid comfort

Equipped with retractable and foldable footrests, desktype armrests that allow easy access to a table and comfortable long armrest padding. Small side wheels on the footplates will protect your walls and furniture while moving around.



Anatomic version (optional) Max. users weight: 150 kg

Colours

C20
blue

Technical data

Total width	560/610/670/720 mm
Total height	920 mm
Total length	1110 mm
Folded width	290/290/310/350 mm
Seat width	390/440/500/550 mm
Height armrests	230 mm

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the model **28** complies with the highest norms, laid down by certified test institutes in most European countries.



Seat height	500 mm
Height backrest	420 mm
Seat depth	430 mm
Weight	+/- 20 kg
Maximum user weight	(150 kg with anatomic cushions) 110 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer