

Eclips X4 90°

90° reclinable backrest



The knowledge and expertise we used to create the Eclips X4 lead to the birth of the Eclips X4 90°.

Position changes throughout the day are beneficial for bloodflow, pressure relief or even just for comfort. With identical features to Eclips X4, Eclips X4 90° offers the additional option of fully reclining the backrest. Add the optional B14 backrest extension and the BZ7 legrests to create maximum comfort and space to lay down. The position of the rear wheel has been shifted to the back to keep the wheelchair in perfect balance and provide safety when the user is laying flat.



VERMEIREN
we care for you

Standard

- Height adjustable handles
- Inclination of the back by hydraulic springs for more comfort during the tilting
- Internal and external collapsible footplates
- Foldable and removable armrests
- Rear wheels with Quickrelease system

Colour



Different adjustment possibilities for the seat height and seat depth



Removable pushbar for a better stability



*Quick release wheels for rapid disassembly
Tipping piece at the rear*



Adjustment of the backrest by means of gas spring for a smooth adjustment

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Eclips X4 90°** complies with the highest norms, laid down by certified test institutes in most European countries.

Technical data

■ Total width	590/620/640/660/680/700/740 mm
■ Total height	910-1015 mm
■ Total length	1000-1210 mm
■ Folded width	300 mm
■ Seat width	390/420/440/460/480/500/550 mm
■ Height armrests	235 mm

■ Seat height	390-520 mm
■ Height backrest	420-460 mm
■ Seat depth	420-500 mm
■ Weight	19,5 kg
■ Maximum user weight	130 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you