

286B

ROLLATOR

| | |
|-----------------------|-------------------|
| Width | 605 mm |
| Length (opened) | 645 mm |
| Seat width | 400 mm |
| Seat height | 610 mm |
| Height handles | 770 - 1010 mm |
| Brakes | multifunctional |
| Weight | 10 kg |
| Incl. basket and tray | 11 kg |
| Max. users weight | 120 kg |
| Max. load basket | 5 kg |
| Wheels | 190x47 |
| Colour | C59 - Carbon grey |

Subject to technical alterations. Measurements in mm, weight in kg.
Deviations in sizes and weight (15 mm and 1,5 kg) are possible.



N.V. Vermeiren N.V. - Vermeirenplein 1-15 - B-2920 Kalmthout - Belgium

Tel.: +32 (0)3 620 20 20 - Fax: +32 (0)3 666 48 94 - info@vermeiren.be - website: www.vermeiren.com

Images for illustration purposes only. Certain features may be optional.

Opening up new horizons in mobility