



Detail A

- . Add loctite nr.243 on nr.34, nr.44 and nr.45
- . Lubricate with grease : nr.11 and nr.19
- . Lubricate with grease : nr.18 and nr.19
- . Use PP02203 KLEJ RED-SABA FOAMSPRAY 2203/2916 to glue nr.52
- . Tighten with torque wrench 25Nm W600135: nr.64 and Nr.27

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

8001500
Plan No.
E00078514/A.1.1

EAGLE 625

Scale:	1:10	Drafted	Verified	Confirmed
Date	28/03/2024	28/03/2024	28/03/2024	28/03/2024
Initials	TT	SD	PJ	



REPLACES E000526

Material			
No.	Plan	Article	Quantity
1	X001016	1001682	5
2	P00074772	2098400	4
3	X004122	1692500	1
4	X004123	1692600	1
5	X004124	1692700	1
6	P00078139	G105000	2
7	X004288	1903587	1
8	X004338	1708700	1
9	X004341	1708800	2
10	X004353	1903588	1
11	X004699	1905496	1
12	P00074779	2089900	2
13	∅12x70	1905630	1
14	Y014084	16897..	1
15	Y014090	16898..	1
16	Y014091	16899..	1
18	P00032386	20824..	1
19	Y016461	17349..	1
20	Y014110	16902..	1
21	P00074771	2098354	1
22	X004720	1905571	1
23	∅25x12.2x0.7	1993500	2
24	∅4x14	1001584	4
25	M8 DIN934	1001671	5
26	M10 DIN985	1001683	4
27	M12 DIN985	1001684	1
28	M8	1006835	4
30	M5x20 ISO7380	1998368	2
31	M10x23	1903086	1
32	M10x80 DIN7991	1905617	2
33	X001398	1009689	1
34	M10x16	1903370	2
35	M12	1003294	1
36	DIN471 17.5-1.2	1904807	1

Material			
No.	Plan	Article	Quantity
37	---	1903402	2
38	---	1903403	2
39	M10x70 DIN7991/8.8	1903406	2
40	---	1918030	1
41	---	1918031	1
42	10 kN	1918028	1
43	---	1903578	1
44	M10x90 DIN7991	1921197	2
45	M8x16	1991858	2
46	M8x35	1993236	2
47	M8x40	1008739	2
48	∅10x50	1998885	1
49	∅25 M10x1	1999935	1
50	---	1918029	1
51	---	1918032	1
52	X004580	1904798	3
53	---	T001052	2
54	---	1904568	2
55	X004612	1904946	1
56	X004632	1905205	1
58	Y016681	F006290	1
60	X004738	1905658	1
61	---	1900641	2
62	∅5x18	1001589	1
64	M12 DIN934	1001673	1
67*	---	1918034	1
68	---	1906233	1
69	M8x25 DIN912/8.8	1001499	1
70	M10 DIN439B	1903216	2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

8001500

Plan No.

E00078514/A1.1

EAGLE 625

Scale: 1:10	Drafted	Verified	Confirmed
Date	28/03/2024	28/03/2024	28/03/2024
Initials	TT	SD	PJ



* Optional charger, not depicted