



Sheet 1/2

2991639

E000818 / K

BACK FRAME L34 TIMIX / FOREST 3 / SIGMA
 RUGFRAME L34 TIMIX / FOREST 3 / SIGMA

Plan No.
 E000818 / K

Scale: 1:4	Drafted	Verified	Confirmed
Date	27/06/2012	15/12/2023	15/12/2023
Initials	SD	PM	MK



Material				Material			
No.	Plan	Article	Quantity	No.	Plan	Article	Quantity
1	---	1908883 (40A)	1	14	X001119	1001554	4
	---	1908877 (40B)	1	15	M8x16 DIN912	1001498	4
	---	1908880 (40C)	1	16	30x20	1997922	2
	---	1906939 (40D)	1	17	M8 DIN125	1003199	2
	---	1915818 (40E)	1	18	ACT Linx	1912307	1
	---	2995303 (40F)	1	19	M4x30 DIN7985	1001487	2
	---	2998162 (40G)	1	20	M4 DIN985	1001678	2
	---	1908884 (45A)	1	21	---	T001002	1
	---	1908878 (45B)	1	22	P00034376	1918563	2
	---	1908881 (45C)	1				
	---	1906941 (45D)	1				
	---	1916611 (45F)	1				
	---	2998163 (45G)	1				
	---	1908885 (50A)	1				
	---	1908879 (50B)	1				
	---	1908882 (50C)	1				
	---	1906943 (50D)	1				
	---	1916417 (50F)	1				
	---	2998164 (50G)	1				
2	X004944	1908268 (40)	1				
	X004983	1908265 (45)	1				
	X004982	1908267 (50)	1				
3	Y017420	1775754 (40)	1				
	Y018139	1769954 (45)	1				
	Y018480	1776454 (50)	1				
4	M6x20 ISO7380	1006999	10				
5	Y017415	1775354 (40)	1				
	Y017991	1769854 (45)	1				
	Y018477	1776154 (50)	1				
6	Y018564	1774200	2				
7	X004991	1773000	2				
8	X001016	1001682	4				
9	---	1908512	2				
10	M12 DIN125	1993362	2				
11	M8x25 DIN7991	1995598	2				
12	M8x20 ISO7380	1992349	2				
13	B58	2993352	1				

- A: Black skai - Blue lining
- B: Black skai - Light green lining
- C: Black skai - Red lining
- D: Grey soft / Black skai
- E: Black skai
- F: Black Cordura
- G: Black soft

Sheet 2/2

2991639

E000818 / K

BACK FRAME L34 TIMIX / FOREST 3 / SIGMA				Plan No.
RUGFRAME L34 TIMIX / FOREST 3 / SIGMA				E000818 / K
Scale: 1:4	Drafted	Verified	Confirmed	
Date	27/06/2012	15/12/2023	15/12/2023	
Initials	SD	PM	MK	