

Left side only depicted

## SIDEFAME ECL X2 KIDS

1/2

Plan No.

EC00450/C

Scale:9:50

Drafted

Verified

Confirmed

Date

19/10/2018

11/11/2020

Init.

Mars

AL




EC00450/C

| Material |              |         |        |
|----------|--------------|---------|--------|
| No.      | Plan         | Article | Number |
| 1        | YC00428(L)   | 1C795   | 1      |
|          | YC00429(R)   | 1C796   | 1      |
| 2        | YC01361      | 1C788   | 2      |
| 3        | YC04035      | 1C97704 | 2      |
| 4        | ---          |         |        |
| 5        | M5 DIN985    | 1001679 | 8      |
| 6        | ---          | 1008492 | 20     |
| 7        | ---          | 1001495 | 6      |
| 8        | ---          | 3991658 | 2      |
| 9        | ---          | 3990684 | 2      |
| 10       | X004689      | 1734600 | 4      |
| 11       | M5X20 DIN912 | 1003907 | 8      |
| 12       | ---          | 1003109 | 4      |
| 13       | ---          | 1001564 | 6      |
| 14       | ---          | 1006850 | 4      |
| 15       | X004309      | 1904133 | 4      |
| 16       | X004818      | 1906184 | 2      |
| 17       | X004415      | 17193   | 2      |
| 18       | ---          | 1903903 | 2      |
| 19       | ---          | 1001610 | 2      |
| 20       | ---          | 1001684 | 2      |
| 21       | ---          | 1992559 | 4      |
| 22       | X000188      | G078500 | 2      |

|    |                 |         |   |
|----|-----------------|---------|---|
| 23 | XC00683         | 3992963 | 2 |
| 24 | XC00014         | 3990031 | 2 |
| 25 | X001016         | 1001682 | 2 |
| 26 | ---             | 3993060 | 2 |
| 27 | ---             | 3990041 | 2 |
| 28 | X004973         | 1907519 | 4 |
| 29 | ---             | 1992835 | 4 |
| 30 | ---             | 1997781 | 2 |
| 31 | XC00037         | 3990106 | 2 |
| 32 | M6x55 DIN912    | 1908726 | 4 |
| 33 | --              | 3990132 | 2 |
| 34 | WHEEL 20" BL PU | 1915557 | 2 |
| 35 | X004469         | 1905045 | 2 |
| 36 | ---             | 3990117 | 2 |
| 37 | ---             | 1001680 | 2 |
| 38 | X005565         | 1912287 | 2 |
| 39 |                 |         |   |
| 40 |                 |         |   |
| 41 | YC02165(L)      | 1D01854 | 1 |
|    | YC02166(R)      | 1D01954 | 1 |
| 42 |                 |         |   |
| 43 |                 |         |   |

EC00450/C

|                       |            |            |           |   |
|-----------------------|------------|------------|-----------|---|
| SIDEFRAME ECL X2 KIDS |            |            |           | Plan No.  |
|                       |            |            |           | EC00450/C   |
| Scale:9:50            | Drafted    | Verified   | Confirmed |  |
| Date                  | 19/10/2018 | 11/11/2020 |           |   |
| Init.                 | Mars       | AL         |           |   |

2/2