


2993391

Sheet 1/2

E001058 / C

FOREST 3+ REAR LOWER FRAME				Plan No.
				E001058 / C
Scale: 1:5	Drafted	Verified	Confirmed	
Date	23/11/2016	01/03/2022	01/03/2022	
Initials	SD	SD	PJ	


Material							
No.	Plan	Article	Quantity	No.	Plan	Article	Quantity
1	X001016	1001682	4	18	M10 DIN985	1001683	2
2	X002282	1992854	4	19	M8x45 DIN912	1003110	4
3	X002792	1555954	4	20	M6x12 ISO7380	1007000	9
4	X003403	1617400	4	21	M5x35 DIN912	1007961	2
5	X004778	1908255	4	22	M12 DIN125	1993362	2
6	X005576	1912735	1	23	---	1906791	2
7	X005609	1912734	1	24	---	1906792	2
8	X005637	1912761	1	25	---	1908165	2
9	X005638	1912762	1	26	M5x10	1909204	2
10	Y020969 DP	20363..	1	27	80A	1909646	2
11	Y021565	2025854	1	28	150A	1910488	1
12	Y022051	2025954	1	29*	---	1912309	1
13	M6 DIN125	1000973	6	30	LINX 120A	1912415	1
15	M5 DIN125	1001550	2	31	M6	1912736	3
16	D28	1001606	2	32	M10x65 DIN912	1912777	2
17	M5 DIN985	1001679	2	33	---	1999843	2
				34	M5x10 ISO7380	2991862	4

* Option SE50

Sheet 2/2

2993391

E001058 / C

FOREST 3+ REAR LOWER FRAME				Plan No.
				E001058 / C
Scale: 1:5	Drafted	Verified	Confirmed	
Date	23/11/2016	01/03/2022	01/03/2022	
Initials	SD	SD	PJ	